
FlowflexTM

SARS-CoV-2

Antigen Rapid Test (Self-Testing)

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The Global Leader of Rapid Test with **26** Years Experience.

Fast

Accurate

Easy to Use

Reliable

CE Marked

Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)

The SARS-CoV-2 Antigen Rapid Test is a lateral flow test for the qualitative detection of the nucleocapsid antigen from SARS-CoV-2 in anterior nasal swab specimens directly from individuals suspected of COVID-19 within the first seven days of the onset of symptoms. The test can also test specimens from individuals without symptoms.

- **Specimen:** Anterior nasal swab specimens
- **Test Time:** Results at 15 min.
- **Shelf life:** 24 months
- **Storage temperature:** 2-30°C
- **Accuracy:** 98.8%
- **Sensitivity:** 97.1%
- **Specificity:** 99.5%

Test Procedure


1. Test Preparation

A. Open your test kit:

B. You should have:

Note: For 1T and 5T, the hole is on the kit box.

Tube Holder
(for 25T only)

Disposable Swab

Test cassette

Extraction Buffer Tube

Package Insert

Waste Bag

C. Wash or sanitize your hands.

D. Read the instructions.


E. Check the expiration date.

F. Open the pouch.


G. Check the cassette


2. Specimen Collection


COLLECTION BY AN ADULT CAREGIVER


SELF COLLECTION

Insert the entire absorbent tip of the swab into one nostril. Using gentle rotation, push the swab less than 2.5 cm from the edge of the nostril.

Rotate the swab 5 times brushing against the inside of the nostril.

Remove the swab and insert it into the other nostril.


Remove swab from the nostril.


3. Sample Preparation


Unscrew the dropper cap.

Insert the swab into the tube and swirl for 30 seconds.

Rotate the swab 5 times while squeezing the side of the tube.


Remove the swab while squeezing the tube.

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BIBLIOGRAFIA

- Shuo Su, Gary Wong, Welfeng Shi, et al. Epidemiology, Genetic recombination, and pathogenesis of coronaviruses. Trends in Microbiology, June 2016, vol. 24, No. 6: 490-502
- Susan R. Weiss, Julian L. Leibowitz, Coronavirus Pathogenesis, Advances in Virus Research, Volume 81: 85-164

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Tamponi monouso

Jiangsu Changfeng Medical Industry Co., Ltd.
Touqiao Town, Guangling District, Yangzhou
225109, Jiangsu, P.R. China

O

Jiangsu HanHeng Medical Technology Co., Ltd.
16-B4, #1 North Qingyang Road, Tianning District,
Changzhou, 213017 Jiangsu P.R. China

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Flowflex™

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

- Contents:
- 25 Test Cassettes
 - 25 Extraction Buffer Tubes
 - 25 Disposable Swabs
 - 25 Waste Bags
 - 1 Tube Holder
 - Package Insert

COVID-19 Test

A rapid test for the detection of SARS-CoV-2 nucleocapsid antigens in anterior nasal swab specimens.
For in vitro diagnostic use only. For self-testing.

Flowflex™

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

- Contents:
- 25 Test Cassettes
 - 25 Extraction Buffer Tubes
 - 25 Disposable Swabs
 - 25 Waste Bags
 - 1 Tube Holder
 - Package Insert

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
YYYY-MM-DD PDV-XXXXXXXXX
LCA9855-01

Attention: By approving the enclosed design draft, you (ACON's Customer) accept all responsibility for the accuracy of the design. If an error is detected following the printing or manufacturing of a material, you (ACON's Customer) are responsible for the cost of any inventory which is deemed unsuitable for sale and/or not compliant with local regulations.

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US	OUS	DOMESTIC	OTHER
Description	Flowflex Bio CE0123 S-2 Ag Test kit box (Self-Testing)-1T	Part Number	1130522601
Printing Contents	✓	L Number	✓
Designer	nafei.chen	Design Date/Version	May 14, 2021 / H
Artwork checked by		Material	350g白卡+水性UV
Approved by Customer		Approved by Marketing/Sales	
Approved by P.M.T.		Approved by QA	Effective Date

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Product Service

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6)
(Devices for self-testing)

No. V9 042074 0032 Rev. 00

Manufacturer:

Acon Biotech (Hangzhou) Co., Ltd.

No.210 Zhenzhong Road
West Lake District
310030 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

Product:

In Vitro diagnostic devices for self testing

The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex III (6). The design of the devices conforms to the requirements of this Directive. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V9_042074_0032_Rev_00

Report No.:

SH2110605

Valid from:

2021-05-14

Valid until:

2024-05-26

Date,

2021-05-14

Christoph Dicks
Head of Certification/Notified Body

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6)
(Devices for self-testing)

No. V9 042074 0032 Rev. 00

Model(s): SARS-CoV-2 Antigen Rapid Test (Self-Testing)

Facility(ies): Acon Biotech (Hangzhou) Co., Ltd.
No.210 Zhenzhong Road, West Lake District, 310030 Hangzhou,
PEOPLE'S REPUBLIC OF CHINA

Model name:

Model number:

Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-11855
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118L5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118T5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118A5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118U5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118V5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118Q5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118M5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118N5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118W5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118P5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118Y5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118Z5
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118R5

SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-11853E
SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118A3E
SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118Q3E
SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118M3E
SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118P3E
SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118R3E

Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-11857
Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118A7
Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118Q7
Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118M7
Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118P7
Quik Check SARS-CoV-2 Antigen Rapid Test (Self-Testing)	REF L031-118R7

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6)
(Devices for self-testing)

No. V9 042074 0032 Rev. 00

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-11853J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118A3J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118Q3J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118M3J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118P3J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118R3J

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-11853D

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118A3D

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118Q3D

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118M3D

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118P3D

SARS-CoV-2 Antigen Rapid Test (Self-Testing)

REF L031-118R3D

Declaration of Conformity

ACON Biotech (Hangzhou) Co., Ltd.
No.210 Zhenzhong Road, West Lake District,
Hangzhou, P.R. China, 310030

**We declare under our sole responsibility that the
in vitro diagnostic device:**

SARS-CoV-2 Antigen Rapid Test (Self-Testing)


**classified as self-testing according to the Annex II of the directive 98/79/EC,
meets all the provisions of the directive 98/79/EC on *in vitro*
diagnostic medical devices which apply to it.**

**This declaration is according to Annex III.6 of the Directive and thus is
based on approval by the notified body
TÜV SÜD Product Service GmbH, Ridlerstraße 65
80339 MÜNCHEN, Germany, notified under
No. 0123 to the EC Commission.**

Authorized Representative:
MedNet GmbH
Borkstrasse 10
48163 Muenster, Germany

This declaration is valid until expiration of EC certificate
No. V9 042074 0032 Rev.00
Expiration Date: 2024-05-26

Signed this 17 day of May, 2021
in Hangzhou, China

Junny You
International Regulatory Affairs Senior Director
ACON Biotech (Hangzhou) Co., Ltd.

ACON BIOTECH (HANGZHOU) CO., LTD.
No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R. China, 310030

FORMATI – CODICI EAN ARTICOLI					
Tests/Kit	Kit Box Dimension	Carton Dimension	Tests/Carton	EAN on Kit Box	EAN on Carton
1T	120*65*27	565*390*365	300	6921756492427	6921756492670
5T	165*80*45	520*395*440	600	6921756492434	6921756492687
25T	210*125*93	445*400*410	600	6921756492533	6921756492694

CODICI ARTICOLO			
Product	REF No.	Count	Tube Type
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-11855	1T	Blue cap tube
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-118A5	5T	Blue cap tube
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-118Q5	25T	Blue cap tube
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-118M5	1T	Foil sealed cap
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-118P5	5T	Foil sealed cap
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-118R5	25T	Foil sealed cap

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FARMADATI ITALIA Srl