VivaDiag™ Pro

SARS-CoV-2 Ag Rapid Test

Enhanced sensitivity

Test Procedure

Specimen Collection

Nasal swab


Oropharyngeal swab


Nasopharyngeal swab

Solution Preparation

1. Prepare the specimen

2. Apply the extracted specimen

3. Read the test result

Detection


Result Interpretation

Positive

Negative

Invalid

Note: The VivaDiag™ Pro SARS-CoV-2 Ag Rapid Test has ONLY been designed to act as a supplementary test for suspected cases of negative coronavirus nucleic acid detection or in conjunction with nucleic acid detection in the diagnosis of suspected cases. Results from nucleocapsid protein antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 (COVID-19) infection or to inform infection status.

VivaChek™

VivaDiag™ Pro SARS-CoV-2 Ag Rapid Test is for the rapid, qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in human nasal swab, oropharyngeal swab or nasopharyngeal swab specimen. The test is for *in vitro* diagnostic use only.

Specification

Test Principle	Immunochromatography
Sample Type	Nasal swab, oropharyngeal swab or nasopharyngeal swab
Sample Volume	3 drops
Test Time	15 min
Operation Temperature	15-30°C
Storage Temperature	2-30°C
Shelf Life (Unopened)	24 months

Not Available for Sales in the US

NOT FOR AT-HOME TESTING
FOR PROFESSIONAL USE ONLY!

VivaChek
VivaChek Biotech (Hangzhou) Co., Ltd.
level 2, Block 2, 146 East Chaofeng Rd.,
Yuhang Economy Development Zone,
Hangzhou, 311100, China
www.vivachek.com

Number:1626001501
Effective date: 2021-09-13
© 2021 VivaChek Biotech (Hangzhou) Co., Ltd.

REF VCD16-01-011/VCD16-01-012/VCD16-01-013/VCD16-01-014/ VCD16-01-015/VCD16-01-016/VCD16-01-017	Italiano
--	----------

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FONTI

1. Coronaviridae Study Group of the International Committee on Taxonomy of Viruses. The species severe acute respiratory syndrome-related coronavirus: classifying 2019-nCoV and naming it SARS-CoV-2 [J]. Nature Microbiology, 5, 536-544 (2020).
2. Perlman, S. Netland, J. Coronaviruses post-SARS: update on replication and pathogenesis. Nature Reviews Microbiology 7, 439-450, doi: 10.1038/nrmicro2147 (2009).
3. Lauer SA, Grantz KH, Bi Q, et al. The Incubation Period of Coronavirus Disease 2019 (COVID-19) From Publicly Reported Confirmed Cases: Estimation and Application. Ann Intern Med. 2020; 172 (9): 577-582. doi: 10.7326/M20-0504.
4. B. Korber et al. Tracking Changes in SARS-CoV-2 Spike: Evidence that D614G Increases Infectivity of the COVID-19 Virus. Cell, vol. 182, no. 4, pp. 812-827.e19, Aug. 2020, doi: 10.1016/j.cell.2020.06.043.
5. Investigation of novel SARS-CoV-2 variant: Variant of Concern 202012/01, GOV.UK.
6. New SARS-CoV-2 variant, GOV.UK.
7. D. C. Dinesh et al. Structural basis of RNA recognition by the SARS-CoV-2 nucleocapsid phosphoprotein. PLoS Pathog. vol. 16, no. 12, p. e1009100, Dec. 2020, doi: 10.1371/journal.ppat.1009100.
8. J. Marién et al. Evaluating SARS-CoV-2 spike and nucleocapsid proteins as targets for antibody detection in severe and mild COVID-19 cases using a Lumines bead-based assay. Journal of Virological Methods, vol. 288, p. 114025, Feb. 2021, doi: 10.1016/j.jviromet.2020.114025.

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Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016 & EN ISO 13485:2016

This is to certify that: VivaChek Biotech (Hangzhou) Co., Ltd.
Level 2, Block 2
146 East Chaofeng Rd.
Yuhang Economy Development Zone
Hangzhou
Zhejiang
311100
China

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Holds Certificate No: **MD 726764**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 & EN ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In Vitro Diagnostic Test Kits (Colloidal Gold).

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For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-07-24

Latest Revision Date: 2020-07-24

Effective Date: 2020-07-24

Expiry Date: 2023-07-23

Page: 1 of 1

...making excellence a habit.™

Evaluation of VivaDiag Pro SARS-CoV-2 Ag Rapid Test

Purpose

The purpose of the study was to evaluate the sensitivity and performance of VivaDiag Pro SARS-CoV-2 Ag Rapid Test by VivaCheck Biotech (Hangzhou) Co., Ltd. with SARS-CoV-2 Nucleocapsid protein and with negative samples.

Materials

VivaDiag Pro SARS-CoV-2 Ag Rapid Test Ref VCD16-01-011, manufacturer VivaCheck Biotech (Hangzhou) Co., Ltd.

Recombinant SARS-CoV-2 Nucleocapsid Protein REF 710081, manufacturer Oy Medix Biochemica AB

Negative nasopharyngeal sample from presumably healthy individual 1 pcs

Method

Recombinant nucleocapsid protein was diluted to the buffer provided by Kit. The interpretation of test results for positive samples between 0,0125 ng/ml and 10 ng/ml were made with IDA solution and visually. All samples were evaluated based on the instruction for use provided by manufacturer. Results were interpreted 15 minutes after adding the sample to each cassette.

Results

Sample	VivaDiag PRO visual result	VivaDiag PRO IDA result
10 ng/ml	++	Positive
1 ng/ml	+	Positive
0,1 ng/ml	+	Positive
0,05 ng/ml	+	Positive
0,0125 ng/ml	+	Positive
0,0076 ng/ml	+/-	Negative
0,00625 ng/ml	+/-	Negative
Buffer	-	Negative
Negative sample 1	-	Negative

Conclusion

VivaDiag Pro SARS-CoV-2 Ag Rapid Test is having exceptionally high analytical sensitivity. Sample having nucleocapsid concentration of 0,0125 ng/ml (12,5 pg/ml) was measured positive both visually and with IDA platform. Negative samples were found to be negative with both visually and with IDA platform.

Based to our evaluation VivaDiag Pro SARS-CoV-2 Rapid Antigen Test is having exceptionally high analytical sensitivity and good specificity.

Anidiagnostics Oy

Karkkila, Finland February 23rd ,2021

Mika Saramäki

CEO

Certificate

Quality Management System EN ISO 13485:2016

Registration No.: SX 2135127-1

Organization: VivaChek Biotech (Hangzhou) Co., Ltd.
Level 2, Block 2, 146
East Chaofeng Rd., Yuhang Economy Development Zone,
Hangzhou,
311100 Zhejiang
P.R. China

Scope: Design and Development, Manufacture and Distribution of Blood Glucose Monitoring Systems (Blood Glucose Meters, Blood Glucose Test Strips and Blood Glucose Control Solutions), In Vitro Diagnostic Test Kits and Analyzers used in the detection of Cardiac Markers, Tumor Markers, Infection, Inflammation, Renal Function, Diabetes, Hormone, Drug of Abuse, Anemia and Vitamins based on fluorescence immunochromatography; Multi-function Monitoring Systems (including Multi-function Meters, Blood Ketone Test Strips, Blood Ketone Control Solutions, Uric Acid Test Strips, Uric Acid Control Solutions, Blood Glucose Test Strips, and Blood Glucose Control Solutions)

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the abovementioned standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.: 244312018-200
Effective date: 2021-03-12
Expiry date: 2024-03-02
Issue date: 2021-03-12

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The scope of certification also covers the following:

No.	Facility	Scope
/01	VivaChek Biotech (Hangzhou) Co., Ltd. Level 2, Block 2, 146 East Chaofeng Rd., Yuhang Economy Development Zone, Hangzhou, 311100 Zhejiang P.R. China	Design and Development, Manufacture of Blood Glucose Test Strips Blood Glucose Control Solutions, In Vitro Diagnostic Test Kits used in the detection of Cardiac Markers, Tumor Markers, Infection, Inflammation, Renal Function, Diabetes, Hormone, Drug of Abuse, Anemia and Vitamins based on fluorescence immunochromatography, Blood Ketone Test Strips, Blood Ketone Control Solutions, Uric Acid Test Strips, Uric Acid Control Solutions

Report No.: 244312018-200
Effective date: 2021-03-12
Expiry date: 2024-03-02
Issue date: 2021-03-12

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
The scope of certification also covers the following:

- /02 VivaChek Biotech (Hangzhou) Co., Ltd.
1/2/3 F, Building 1,
16 East Zhenxing Rd., Yuhang Economy
Development Zone,
Hangzhou,
311100 Zhejiang
P.R. China
- Design and Development, Manufacture of
Blood Glucose Meters, Multi-function Meters,
In Vitro Diagnostic Analyzers used in the
detection of Cardiac Markers, Tumor
Markers, Infection, Inflammation, Renal
Function, Diabetes, Hormone, Drug of
Abuse, Anemia and Vitamins based on
fluorescence immunochromatography
Distribution of Blood Glucose Monitoring
Systems (Blood Glucose Meters, Blood
Glucose Test Strips and Blood Glucose
Control Solutions); In Vitro Diagnostic Test
Kits and Analyzers used in the detection of
Cardiac Markers, Tumor Markers, Infection,
Inflammation, Renal Function, Diabetes,
Hormone, Drug of Abuse, Anemia and
Vitamins based on fluorescence
immunochromatography; Multi-function
Monitoring Systems (including Multi-function
Meters, Blood Ketone Test Strips, Blood
Ketone Control Solutions, Uric Acid Test
Strips, Uric Acid Control Solutions, Blood
Glucose Test Strips, and Blood Glucose
Control Solutions)

Report No.: 244312018-200
Effective date: 2021-03-12
Expiry date: 2024-03-02
Issue date: 2021-03-12

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Product Name and Model

VivaDiag™ Pro SARS-CoV-2 Ag Rapid Test VCD16
(The product has been registered and sold in Europe)

Classification:

Other device, not in annex II, not for self-testing, not for performance evaluation.

Conformity assessment procedure: ANNEX III, 98/79/EC

We hereby declare that the above mentioned products meet the COUNCIL DIRECTIVE 98/79/EC and applicable standards. All supporting documentations are retained in the manufacturer and EU representative.

General applicable standards:

DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices.

Hangzhou, China, 5th Feb, 2021

Place, Date of issue

Elaine Huang
Regulatory Affairs Manager

VivaChek Biotech (Hangzhou) Co., Ltd

Level 2, Block 2, 146 East Chaofeng Rd, Yuhang Economy Development Zone, Hangzhou, Zhejiang 311100, China
Tel: +86-571-89182700 Fax: +86-571-89182733 Email: info@vivachek.com www.vivachek.com

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Lotus NL B.V.