ALL TEST

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【BIBLIOGRAFIA】

1. Steier JA, P Bergsjö, OL Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy, *Obstet Gynecol.* 1984; 64(3): 391-394
2. Dawood MY, BB Saxena, R Landesman Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma, *Obstet. Gynecol.* 1977; 50(2): 172-181
3. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production of human chorionic gonadotropin by neoplasms", *Ann. Intern Med.* 1973; 78(1): 39-45

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Manufacturer:

Name: HANGZHOU ALLTEST BIOTECH CO., LTD.

Address: #550, Yin Hai Street, Hangzhou Economic & Technological Development Area, Hangzhou -310018, P.R. China

European Representative:

Name: MedNet EC-REP GmbH

Address: Borkstrasse 10, 48163 Muenster, Germany

Product Name: HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)

Cat. No.: See Attachment 1

Analyte: Human Chorionic Gonadotropin in human Urine

Model: See Attachment 1

Classification: Self-testing of IVDD 98/79/EC

Conformity Assessment Route: IVDD 98/79/EC Annex IV

EDMA Code: 12 70 05 02 00

We, HANGZHOU ALLTEST BIOTECH CO., LTD., herewith declare that we are exclusively responsible for this declaration of conformity. We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.

DIRECTIVES

General applicable directives:

DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices

Standard Applied: EN ISO 13485:2016, EN ISO 14971:2012, EN 13975:2003, EN ISO 18113-1:2011, EN ISO 18113-4:2011, EN 13612:2002/AC: 2002, EN ISO 17511:2003, EN ISO 23640:2015, EN 13641:2002, EN 13532:2002, EN ISO 15223-1:2016


Notified body: TUV SUD Product service GmbH, Ridlerstrasse 65, 80339 Munich, Germany (0123)

(EC) Certificate(s): V1 095123 0008 Rev.04

Expire date of the Certificate: 2025-05-26

Place, Date of First issue of DOC: in Hangzhou on 2016-09-14

The Date of Issue of DOC on 2022-04-05

Signature: 

Name: GAO FEI (Position: General Manager)

Hangzhou AllTest Biotech Co.,Ltd.
#550, Yin Hai Street,
Hangzhou Economic & Technological
Development Area,
Hangzhou -310018, P.R. China

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Attachment 1

Catalog NO.	Product Name	Analyte	EDMA Code	Model
FHC-101H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Dipstick
FHC-102H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Cassette
FHC-103H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Midstream
FHC-F101H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Dipstick
FHC-F102H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Cassette
FHC-F103H	HCG (Human Chorionic Gonadotropin) Pregnancy Rapid Test (Urine)	Human Chorionic Gonadotropin in human Urine	12 70 05 02 00	Midstream

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Certificate

No. Q5 095123 0007 Rev. 04

Holder of Certificate: **Hangzhou AllTest Biotech Co., Ltd.**
550#, Yin Hai Street
Hangzhou Economic and Technological Development Area
310018 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: **Design and Development, Production and Distribution of In Vitro Diagnostic Reagents, Control Material and Instruments for Clinical Chemistry, Immunochemistry (Immunology), Haemostasis, and Infectious Diseases, including Professional Laboratory Use, Near Patient and Self Testing**

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: [www.tuvsud.com/ps-cert?q=cert:Q5 095123 0007 Rev. 04](http://www.tuvsud.com/ps-cert?q=cert:Q5_095123_0007_Rev.04)

Report No.: SH23106403, SH23106403_CN

Valid from: 2023-09-25
Valid until: 2026-09-24

Date, 2023-09-20

Christoph Dicks
Head of Certification/Notified Body

Certificate

No. Q5 095123 0007 Rev. 04

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2016)
DIN EN ISO 13485:2016

Facility(ies): **Hangzhou AllTest Biotech Co., Ltd.**
550#, Yin Hai Street, Hangzhou Economic and Technological
Development Area, 310018 Hangzhou, PEOPLE'S REPUBLIC OF
CHINA

Design and Development, Production and Distribution of In Vitro
Diagnostic Reagents, Control Material and Instruments for Clinical
Chemistry, Immunochemistry (Immunology), Haemostasis, and
Infectious Diseases, including Professional Laboratory Use, Near
Patient and Self Testing

Hangzhou AllTest Biotech Co., Ltd.
#383, Qiaoxin Road, Xiasha Street, Qiantang District, 310018
Hangzhou, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

Design and Development, Production and Distribution of In Vitro
Diagnostic Reagents, Control Material and Instruments for Clinical
Chemistry, Immunochemistry (Immunology), Haemostasis, and
Infectious Diseases, including Professional Laboratory Use, Near
Patient and Self Testing